



Model:
CL-14

Typology:
ILLUMINATOR

Description:
LED ring light illuminator.

Cromaticity coordinates (x,y): 0.320, 0.320 (pure white), 6,300 K.

Optical power: Luminous flux: 240 lumen.
 Illuminance (at 10 cm distance): 8,000 lux (on 3 cm diameter).

Power supply: Brightness adjustment knob. Main power switch 0/1.
 External power supply provided.
 Input: 110/240 Vac 50/60 Hz; output: 12 Vdc 1,500 mA.

Suitable models: LAB, SZM and SZN.

Accessories: Fixing ring for stereo heads.