



Model:

CL-41

Typology:

ILLUMINATOR

Description:

Double flexible arm (L = 60 cm) with double X-LED³.

Chromaticity coordinates (x,y):

0.320, 0.320 (pure white).

Electrical power

3.6 W (each arm).

Optical power:

Luminous flux: 400 lumen.

Illuminance (at 10 cm distance): 170,000 lux.

Aperture angle

25°.

Power supply:

Brightness control potentiometer.

External power supply provided.

Input: 110/240 Vac 50/60 Hz; output: 12 Vdc 1500 mA.