



Coaxial coarse & fine focusing system

Focusing System	Type	Coaxial coarse & fine
	Coarse total travel (mm)	50
	Fine total travel (per single rotation) (mm)	0.2
	Fine graduations	100
	Fine resolution (μm)	2
	Adjustable tension	Yes
	Head holder internal diameter (mm)	76
Accessories Included	User Manual	Digital version (downloadable)

Accessories included	User Manual	Digital version (downloadable)
Product Dimensions	Total height (mm)	100
	Total width (mm)	185
	Total depth (mm)	250
Product Weight	(kg)	1.1